Zeeman Effect

Henry Su, Wen Hong Lou NTU Physics Department

April 4, 2025

Abstract

This report investigates the Zeeman Effect, a phenomenon where spectral lines are split into multiple components in the presence of a magnetic field. The experiment aims to measure the splitting and verify the theoretical predictions.

Contents

| 1 | Introduction | 2 | | |
|---|---|---|--|--|
| 2 | Theory | 2 | | |
| | 2.1 Zeeman Effect | 2 | | |
| | 2.2 Selection Rules | 2 | | |
| | 2.3 Normal and Anomalous Zeeman Effect | 3 | | |
| | 2.4 Fabry-Perot Etalon: Frequency Selection | 4 | | |
| 3 | Experiment Setup and Procedure | | | |
| | 3.1 Experimental Setup | 6 | | |
| | 3.2 Procedure | 7 | | |
| | 3.3 Transversal Observation | 7 | | |
| | 3.4 Calibration Remark | 7 | | |
| 4 | Results | | | |
| 5 | Discussion | | | |
| 6 | Conclusion | | | |

1 Introduction

The Zeeman effect, discovered by Pieter Zeeman in 1896, is the splitting of atomic spectral lines in the presence of a magnetic field. This experiment investigates both the normal and anomalous Zeeman effects, demonstrating how magnetic fields influence atomic energy levels. Using a cadmium lamp and a Fabry-Perot interferometer, we will observe the shifts in wavelength caused by the interaction of electron magnetic moments with an external magnetic field. Additionally, we will calculate the Bohr magneton, a fundamental physical constant that quantifies the magnetic moment of an electron due to its orbital or spin motion. This study enhances our understanding of atomic structure and the fundamental principles of quantum mechanics.

2 Theory

2.1 Zeeman Effect

Zeeman effect is a phenomenon that an atomic spectral line is split into several components in the presence of a magnetic field. The splitting occurs due to the interaction between the magnetic moment of the atom and the external magnetic field. The effect can be classified into two types: normal and anomalous Zeeman effects. We'll get into the details of these two types in the latter section. But for now let us focus on the basic theory of the Zeeman effect. In the presence of weak magnetic field B the perturbed Hamiltonian can be written as follows:

$$H' = -\mu \cdot \mathbf{B} = -\mu_B \mathbf{J} \cdot \mathbf{B}, \quad \mu_B = \frac{e\hbar}{2m_e} = 9.274 \times 10^{-24} \text{J/T}$$
 (1)

where $\vec{\mu}$ is the magnetic moment, μ_B is the Bohr magneton, and **J** is the total angular momentum of the atom. In particular, $\mathbf{J} = \mathbf{L} + 2\mathbf{S}$ where the factor 2 is the g-factor for the electron spin. In the latter on experiments our goal is to calculate Bohr magneton through energy level spliting.

Now let's consider the energy of the perturbed Hamiltonian by firt order perturbation theory. The "good state" of the perturbion is the normalized eigenbasis of the unperturbed Hamiltonian(H_0), and it's characterized mainly by magnetic quantum number m_j and orbital quantum number ℓ . With $\mathbf{B} = \mathrm{B}\hat{z}$, the energy of the perturbed Hamiltonian can be written as:

$$E_{mj} = \langle n\ell j m_j \mid H' \mid n\ell j m_j \rangle \tag{2}$$

$$= \mu_B B \left[1 + \frac{j(j+1) + s(s+1) - l(l+1)}{2j(j+1)} \right] m_j = -g_j \mu_B B m_j.$$
 (3)

where g_i is the Landé g-factor.

2.2 Selection Rules

Dispite knowing physically all the possible energy level existed, naturally not all transitions between state is allowed. We know the selection rule coonstrains the possible transitions between these energy states. There are two types of transitions allowed:

$$\pi$$
-line: $\Delta m_i = \pm 1$ (4)

$$\sigma_{\pm}$$
-line: $\Delta m_i = 0$ (5)

In our experimentm, we detected lights propagating parellel and perpendicular to the magnetic field B. Furthermore, the corrsponding two types of transitions to the light propagatind directions can be summurized in Fig.1




Figure 1: Illustration of the allowed transitions for the Zeeman effect, showing π -lines $(\Delta m_j = 0)$ and σ -lines $(\Delta m_j = \pm 1)$ for light propagating perpendicular to the magnetic field. Note that there is no π -line for light propagating parallel to the magnetic field because the magnetic dipole is oriented along the field. [1]

2.3 Normal and Anomalous Zeeman Effect

The primary distinction between the normal and anomalous Zeeman effects lies in the electron spin. For the normal Zeeman effect, the electron spin s=0, whereas for the anomalous Zeeman effect, $s \neq 0$.

In terms of transitions, the initial state can be represented as $|i\rangle = |\ell_i, s_i, j_i, m_{j_i}\rangle$, and the final state as $|f\rangle = |\ell_f, s_f, j_f, m_{j_f}\rangle$. In our experiment, we use cadmium, and the specific transitions for the normal and anomalous Zeeman effects are as follows:

Normal:
$$3^{1}D_{2} \to 2^{1}P_{1}$$
 or $|3; 2, 0, 2, m_{j}\rangle \to |2; 1, 0, 1, m_{j}\rangle$, (6)

Anomalous:
$$2^{3}S_{1} \to 2^{3}P_{2}$$
 or $|2; 0, 1, 1, m_{j}\rangle \to |2; 1, 1, 2, m_{j}\rangle$. (7)


Indeed this form of transition representation loses the information of quantum number m_j . However, that's a not a big deal since later we only care about the Δm_j . So going back to the enegry difference between transions, we can use eq.3 to calculate the energy difference between the two transitions as:


$$\Delta E = \Delta E_0 + \mu_B B \left[g_{j_f} m_{j_f} - g_{j_i} m_{j_i} \right], \begin{cases} \Delta m_j = m_{j_f} - m_{j_i} = 0, \text{ for } \pi\text{-line} \\ \Delta m_j = m_{j_f} - m_{j_i} = \pm 1, \text{ for } \sigma_{\pm}\text{-line} \end{cases}$$
(8)

with the plotting corresponding quantum number in Eq. 7, we have,

$$\begin{cases}
Normal: g_{j_f} = g_{j_i} = 1, \ \Delta m_j = m_{j_f} - m_{j_i}, \\
Anomalous: g_{j_i} = 2, \ g_{j_f} = \frac{3}{2}
\end{cases}$$
(9)

Then we can summarize the transitions in Fig. 2.





(a) Normal Zeeman Transitions

(b) Anomalous Zeeman Transitions

Figure 2: Transitions we use for cadmium in our experiment. Theoretically, there are nine possible transition states for both effects. However, after passing through the Fabry-Perot Étalon, some transitions fail to exist.[1]

2.4 Fabry-Perot Etalon: Frequency Selection

The Fabry-Perot etalon is a critical component in our experiment, used to select specific frequencies from the light emitted by the cadmium lamp. Since the light source is not monochromatic, the etalon ensures that only the desired spectral lines are transmitted for analysis.

The etalon consists of two parallel, partially reflecting mirrors separated by a distance d. Light waves reflect back and forth between the mirrors, creating an interference pattern. The condition for constructive interference is given by:


$$2\mu t \cos \theta = n\lambda,\tag{10}$$


where t is the separation between the mirrors, θ is the angle of incidence, n is the interference order, μ is the reflective index in quartz, and λ is the wavelength of the transmitted light.


Figure 3a illustrates the working principle of the Fabry-Perot etalon. Light entering the etalon at an angle θ undergoes multiple reflections, and only wavelengths satisfying the interference condition (eq.10) emerge constructively. Thus we have:

$$n = \frac{2\mu t}{\lambda} \cos \theta_n = n_0 \cos \theta_n$$
, where $n_0 = \frac{2\mu t}{\lambda}$. (11)

$$\theta_n = \sqrt{\frac{2(n_0 - n)}{n_0}}\tag{12}$$







(b) Ring radius on the analyser is $r = f\theta$, where f is the focal length of the lens.

And since $n_1 < n_0$ and $n_1 = n_0 \cos \theta$, we can let:

$$n_1 = n_0 - \epsilon \quad ; 0 < \epsilon < 1 \tag{13}$$

And for the pth order of interference, we have:

$$n_p = (n_0 - \epsilon) - (n_p - 1) \tag{14}$$

We can see in figure 3b that the light rays are focused on the analyser. With small angle approximation, we can write the radius of the ring on the analyser as:

$$r = f \tan \theta \approx f\theta \tag{15}$$

Therefore, combining eq.10 and eq.15, we can get the difference on radius of the ring on the analyser as:

$$r_{p+1}^2 - r_p^2 = 2\frac{f^2}{n_0}, \quad \text{with } r_p = \sqrt{\frac{2f^2}{n_0}} \cdot \sqrt{(p-1) + \epsilon}$$
 (16)

Now if there are two different wavelengths λ_a and λ_b in the light source, we can write the n_0 for both wavelengths respect to its corresponding wave number k_a and k_b as:

$$\begin{cases} \epsilon_a &= \frac{2\mu \cdot t}{\lambda_a} = -n_{1,a} = 2\mu \cdot t \cdot k_a - n_{1,a}, \\ \epsilon_b &= \frac{2\mu \cdot t}{\lambda_b} = -n_{1,b} = 2\mu \cdot t \cdot k_b - n_{1,b}. \end{cases}$$
(17)

Then combining eq.16 and eq.17, we can get the difference of wave number k as:

$$\Delta k = k_a - k_b = \frac{1}{2\mu \cdot t} \left[\frac{r_{p+1,a}^2}{r_{p+1,b}^2 - r_{p,a}^2} - \frac{r_{p+1,b}^2}{r_{p+1,b}^2 - r_{p,b}^2} \right]$$
(18)

From the relation in eq.15, we can assume for $n_{0,b} \approx n_{0,a}$ we have:

$$\Delta_a^{p+1,p} = r_{p+1}^2 - r_{p,a}^2 = 2\frac{f^2}{n_{0,a}} = \Delta_b^{p+1,p} = r_{p+1}^2 - r_{p,b}^2 = 2\frac{f^2}{n_{0,b}}$$
(19)

and for all value of p we have:

$$\delta_a^p = r_{p+1,a}^2 - r_{p+1,a}^2 \tag{20}$$

Thus, we rewrite eq.18 as(assuming $\Delta E_0 = 0$):

$$\Delta k = \frac{1}{2\mu \cdot t} \cdot \frac{\delta}{\Delta} \tag{21}$$

From the above we can require the energy difference in terms of Δk :

$$\Delta E = \mu_B B \left[g_{j_f} m_{j_f} - g_{j_i} m_{j_i} \right] = \hbar c \Delta k = \frac{\hbar c}{2\mu \cdot t} \cdot \frac{\delta}{\Delta}$$
 (22)

To summarize, we can see the possible rings in given p in table.1

| | (a) $\Delta(g_j m_j)$ | (b) Possible rings for a fixed order |
|-----------|--|--|
| Normal | $\sigma(\Delta m_j = +1) : \{+1, +1, +1\}$ $\sigma(\Delta m_j = -1) : \{-1, -1, -1\}$ $\pi(\Delta m_j = 0) : \{0, 0, 0\}$ | Transverse: $1 + 2$ rings Longitudinal: $0 + 2$ rings |
| Anomalous | $ \begin{aligned} \sigma(\Delta m_j &= +1) : \{2, \frac{3}{2}, 1\} \\ \sigma(\Delta m_j &= -1) : \{-\frac{3}{2}, -1, -\frac{1}{2}\} \\ \pi(\Delta m_j &= 0) : \{\frac{1}{2}, 0, -\frac{1}{2}\} \end{aligned} $ | Transverse: $3 + 6$ rings Longitudinal: $0 + 6$ rings |

Table 1: (a) The value of $\Delta(g_j m_j)$ for normal and anomalous Zeeman effects, and (b) the corresponding number of rings in the Fabry–Perot interference pattern for a fixed order with the left number representing π -line and right number representing σ_{\pm} .

3 Experiment Setup and Procedure

3.1 Experimental Setup

The electromagnet is placed on a rotating table designed for heavy loads and is mounted with two pole-pieces that have holes, leaving a gap (9–11 mm) large enough for the Cd-lamp. The pole-pieces must be tightened securely. The Cd-lamp is inserted into the gap without touching the pole-pieces and is connected to the power supply for spectral lamps.

The scheme of the set up can be shown in figure.4.




Figure 4: Experimental setup for observing the Zeeman effect. The optical bench is equipped with a lens system, a Fabry-Perot étalon, and a CCD camera for capturing the interference pattern.

3.2 Procedure

- 1. First we open the cadmium light and put on the red light filter make sure that the fringes on the CCD is center at the center of the screen.
- 2. Then we adjust the Fabry-Perot etalon to get the best interference pattern. The etalon is adjusted by rotating the two screws on the side of the etalon. The best interference pattern is achieved when the rings are sharp and well-defined. We also change focals and exposure of the CCD to get the best resolution
- 3. After that we took picture of the polarization entering the CCD at 0° , 90° , and 180°
- 4. Then we rotate the electromagnet by 90°, insert quarter waveplate and replace the red light filter with green light filter for longitudinal observation.
 - 5. We repeat the same procedure as above to get the best interference pattern.
 - 6. After that we took picture of the polarization entering the CCD at 0°, 90°, and 180°

3.3 Transversal Observation

After completing the initial adjustments, the electromagnet is rotated by 90°, and the iris diaphragm is inserted for transversal observation.

3.4 Calibration Remark

For later evaluations, a calibration curve of the magnetic flux density versus the coil current must be recorded. This can be performed using a teslameter. If unavailable, the results shown in Fig. 3 can be used. Note that the flux density measured in the center of the gap (in the absence of the Cd-lamp) should be increased by 3.5% to account for the non-uniform flux distribution within the gap.

4 Results

Present the observed spectral line splitting and compare it with theoretical predictions. Include tables and graphs where necessary.

5 Discussion

Analyze the results, discuss sources of error, and evaluate the agreement between experimental and theoretical values.

6 Conclusion

Summarize the findings and their implications for understanding the Zeeman Effect.

References

References

 $[1] \quad \text{PHYWE Systeme GmbH \& Co. KG. } \textit{Normal and Anomalous Zeeman Effect. 2010}.$